



Characteristics

Capacity	41 lt (1,45 cu ft)
External Dimensions (W x D x H)	cm 60 x 37 x 40
	<i>Internal Dimensions (W x D x H) cm 37x 28,5 x 30 (+ additional compartment 15 x 28,5 x 12,5); Weight kg 17</i>
Structure and Insulation	Atoxic plastic material int./est., thickness 40 mm
Door	Nr 1 blind lid
Key Lock	
Internal equipment	Nr 1 open wire basket dim. cm. 28 x 25 x 25 (W x D x H) Nr.2 cables: for 12/24V and for 110/240V
Controller	Digital electronic thermostat
Alarms	Acoustic and Visual for temperature deviations
Refrigeration Type	Static
Defrost	manual

Technical data

Details	<ul style="list-style-type: none"> • Refrigerant type: R600a CFC free • Mains Voltage 12V/24V; 110V/240V/50Hz-60Hz • Nominal consumption 12 Vcc: 2,58 A • Nominal power consumption: 31 W • Ambient Temperature: +10/+32 °C; humidity 60% • Connectors: 2 Sockets (1 for 12/24V and 1 for 230V)
In accordance with	2014/30/CE, 2014/35/CE and following modifications
Package	Strong Cardboard on wooden pallet
	<i>Dim. cm.64x42x56; Kg.23</i>
	<ul style="list-style-type: none"> • Wooden crate on wooden pallet available: dim. cm.75x50x62; Kg.36

Accessories

Transformer 230V	+10 / -10 °C
-------------------------	--------------

Accessories

Graphic Chart Recorder

Datalogger

FIOCCHETTI SCIENTIFIC S.R.L.

Via Panagulis, 48 - 42045 Luzzara (RE) Italy
Tel. +39 0522 976232- Fax +39 0522 976028
info@fiocchetti.it - www.fiocchetti.it
P.I.00934960352
Registro A.E.E. IT8020000002833