



Characteristics

Capacity	65 lt (2,30 cu ft)
External Dimensions (W x D x H)	cm 70 x 44 x 47
	<i>Internal Dimensions (W x D x H) cm 47 x 37,5 x 35 (+ additional compartment 15 x 37,5 x 17); Weight kg 24</i>
Structure and Insulation	Atoxic plastic material int./est., thickness 40 mm
Door	Nr 1 blind lid
Key Lock	
Internal equipment	Nr 1 open wire basket dim. cm. 46 x 36 x 32 (W x D x H) Nr.2 cables: for 12/24V and for 110/240V
Controller	Digital electronic thermostat
Alarms	Acoustic and Visual for temperature deviations
Refrigeration Type	Static
Defrost	manual

Technical data

Details	<ul style="list-style-type: none"> • Refrigerant type: R600a CFC free • Mains Voltage 12V/24V; 110V/240V/50Hz-60Hz • Nominal consumption 12 Vcc: 3,68 A • Nominal power consumption: 45 W • Ambient Temperature: +10/+32 °C; humidity 60% • Connectors: 2 Sockets (1 for 12/24V and 1 for 230V)
In accordance with	2014/30/CE, 2014/35/CE and following modifications
Package	Strong Cardboard on wooden pallet
	<i>Dim. cm.80x60x63; Kg.32</i>
	<ul style="list-style-type: none"> • <i>Wooden crate on wooden pallet available: dim. cm.85x57x69; Kg.45</i>

Accessories

Transformer 230V	+10 / -20 °C
-------------------------	--------------

Accessories

Graphic Chart Recorder

Datalogger

FIOCCHETTI SCIENTIFIC S.R.L.

Via Panagulis, 48 - 42045 Luzzara (RE) Italy
Tel. +39 0522 976232- Fax +39 0522 976028
info@fiocchetti.it - www.fiocchetti.it
P.I.00934960352
Registro A.E.E. IT8020000002833