

## AF1124W



### Characteristics

<b>External Dimensions (W x D x H)</b>	cm 95 x 65,5 x 91,5
<b>Structure and Insulation</b>	stainless steel inside and outside
<b>Height adjustable Feet</b>	
<b>Kind of ice</b>	flakes
<b>Motor cooling</b>	water
<b>Production Kg/24h (up to)</b>	120
<b>Bin capacity (kg)</b>	40
<b>Weight</b>	kg 78

### Technical data

#### Details

- Refrigerant type: R134A CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption: 12,7 kW/h
- Water consumption l/h: 29,8
- Electrical consumption: 500 W
- Fuse (A): 10
- Cooling Capacity Btuh/W: 1420
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 64 dB
- Ambient Temperature: +18/+43 °C; humidity 60%

<b>In accordance with</b>	2014/30/CE, 2014/35/CE and following modifications
<b>Package</b>	Strong Cardboard on wooden pallet  <i>Dim. cm. 102 x 66 x 107; Kg.84</i> <ul style="list-style-type: none"><li>• <i>Wooden crate on wooden pallet available: dim. cm. 112 x 74 x 115; Kg. 106</i></li></ul>
<b>note</b>	Available also air cooled

### Accessories

Kit for Demineralized Water

60 Hz Frequency

115 V

Wooden crate Package

#### **FIOCCHETTI SCIENTIFIC S.R.L.**

Via Panagulis, 48 - 42045 Luzzara (RE) Italy  
Tel. +39 0522 976232- Fax +39 0522 976028  
info@fiocchetti.it - www.fiocchetti.it  
P.I.00934960352  
Registro A.E.E. IT8020000002833