

## B30A



### Characteristics

**External Dimensions (W x D x H)** cm 39 x 51,5 x 64

**Structure and Insulation** stainless steel inside and outside

**Kind of ice** cubers

**Motor cooling** air

**Production Kg/24h (up to)** 28

**Bin capacity (kg)** 8,5

**Weight** kg 36

### Technical data

#### Details

- Refrigerant type: R134A CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption Kwh/24H: 7 kW/h
- Water consumption l/h: 4
- Fuse (A): 10
- Cooling Capacity Btu/h/W: 2286
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 58 dB
- Ambient Temperature: +10/+32 °C; humidity 60%

**In accordance with** 2014/30/CE, 2014/35/CE and following modifications

**Package** Strong Cardboard on wooden pallet

*Dim. cm.53 x 66 x 89; Kg.44*

- *Wooden crate on wooden pallet available: dim. cm. 59 x 79 x 99; Kg. 68*

**note** Available also water cooled

### Accessories

Kit 4 feet (+135mm H)

60 Hz Frequency

Wooden crate Package

#### **FIOCCHETTI SCIENTIFIC S.R.L.**

Via Panagulis, 48 - 42045 Luzzara (RE) Italy  
Tel. +39 0522 976232- Fax +39 0522 976028  
info@fiocchetti.it - www.fiocchetti.it  
P.I.00934960352  
Registro A.E.E. IT8020000002833