

B40A



Characteristics

External Dimensions (W x D x H) cm 46,7x 57 x 69

Structure and Insulation stainless steel inside and outside

Height adjustable Feet

Kind of ice cubers

Motor cooling air

Production Kg/24h (up to) 40

Bin capacity (kg) 15

Weight kg 35

Technical data

Details

- Refrigerant type: R134A CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption Kwh/24H: 9 kW/h
- Water consumption l/h: 5
- Fuse (A): 10
- Cooling Capacity Btu/h/W: 4094
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 59 dB
- Ambient Temperature: +10/+32 °C; humidity 60%

In accordance with 2014/30/CE, 2014/35/CE and following modifications

Package

Strong Cardboard on wooden pallet

Dim. cm.70 x 66 x 105; Kg.55

- *Wooden crate on wooden pallet available: dim. cm.74 x 74 x 106; Kg.71*

note

Available also water cooled

Accessories

Kit 4 feet (+135mm H)

60 Hz Frequency

Wooden crate Package

FIOCCHETTI SCIENTIFIC S.R.L.

Via Panagulis, 48 - 42045 Luzzara (RE) Italy
Tel. +39 0522 976232- Fax +39 0522 976028
info@fiocchetti.it - www.fiocchetti.it
P.I.00934960352
Registro A.E.E. IT8020000002833