

## B90W



### Characteristics

|  |                                    |
|--|------------------------------------|
| <b>External Dimensions (W x D x H)</b> | cm 70 x 60 x 97 / 109              |
| <b>Structure and Insulation</b>        | stainless steel inside and outside |
| <b>Height adjustable Feet</b>          |                                    |
| <b>Kind of ice</b>                     | cube                               |
| <b>Motor cooling</b>                   | water                              |
| <b>Production Kg/24h (up to)</b>       | 96                                 |
| <b>Bin capacity (kg)</b>               | 50                                 |
| <b>Weight</b>                          | kg 58                              |

### Technical data

#### Details

- Refrigerant type: R134A CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption kWh/24H: 16
- Water consumption l/h: 52
- Fuse (A): 16
- Cooling Capacity Btu/h/W: 2270
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 58 dB
- Ambient Temperature: +18/+43 °C; humidity 60%

|                           |  |
|---------------------------|--|
| <b>In accordance with</b> | 2014/30/CE, 2014/35/CE and following modifications   |
| <b>Package</b>            | Strong Cardboard on wooden pallet<br><br><i>Dim. cm.76 x 63 x 106; Kg.68</i><br><br>• Wooden crate on wooden pallet available: dim. cm.86 x 68 x 111; Kg. 77 |
| <b>note</b>               | Available also air cooled  |

### Accessories

60 Hz Frequency

Wooden crate Package

#### **FIOCCHETTI SCIENTIFIC S.R.L.**

Via Panagulis, 48 - 42045 Luzzara (RE) Italy  
Tel. +39 0522 976232- Fax +39 0522 976028  
info@fiocchetti.it - www.fiocchetti.it  
P.I.00934960352  
Registro A.E.E. IT8020000002833