

BODIES TRANSPORT TROLLEY

Cod. CRR002



Disclaimer:
Images are purely
indicative. The features
and the accessories of the
products may differ
according to model and
application.

General features

External Dimensions (WxDxH)	cm 206x61x81,6
Rollers	4 castors (200 mm diameter), 2 with brakes
Structure	It is made entirely of AISI 304 stainless steel and consists of an electro-welded tubular support structure with no screws or rivets, shaped with a central inclination to contain any liquid.
Maximum load capacity	150 kg

Technical data

Details	The stretcher can be completely extracted thanks to the handles located on head and foot sides.
----------------	---

Accessories

Cover for hydraulic trolley

Wooden crate Package

