



Disclaimer:
Images are purely
indicative. The features
and the accessories of the
products may differ
according to model and
application.

General features

External Dimensions (WxDxH)	cm 206x69x38,5/215
Rollers	2 swivel casters and 2 free ones
Structure	It is constructed in AISI 304 stainless steel tubing and profiles with a double scissors lifting frame. The hydraulic lifting mechanism is powered by a hydraulic pump and is operated via a foot pedal. A flow regulator prevents a too quick lowering of load after lifting, possibly harmful for the operator.
Maximum load capacity	220 kg

Technical data

Details	Resting surface provided with rollers to support the stretcher (not provided) to be inserted or extracted from the mortuary chamber. The trolley is equipped with stops to secure the stretcher, an ergonomic handle and front anti-collision bumpers.
Notes	The equipment is fully compliant with the specific requirements of the norms: UNI EN ISO 12100/1 – 12100/2, UNI EN ISO 954-1, UNI EN ISO 349:2008 including subsequent modifications and integrations.

Accessories

Cover for hydraulic trolley
Wooden crate Package