

General features



| | |
|------------------------------------|---|
| External Dimensions (WxDxH) | cm.215x75x26 / 161 |
| Rollers | 4 steerable solid rubber wheels, 2 of them lockable, 1 enabling locking as regards the direction, d=125 mm with special roller bearings |
| Structure | Precoated steel, with 5 carrying rollers of corrosion-resistant material and plastic coating |
| Weight | 100 kg |

Disclaimer:
Images are purely
indicative. The features
and the accessories of the
products may differ
according to model and
application.

Technical data

| | |
|----------------|--|
| Details | <ul style="list-style-type: none"> Hydraulic unit for foot-actuated vertical adjustment Holding device to secure the body tray Removable slide handle and anti-collision buffer at the front end Equipped with n. 2 U-profiles to accomodate the body tray |
| Note | <ul style="list-style-type: none"> Carrying capacity: max. 250 Kg Suitable for side-lading mortuaries with storage up to 3 levels |

Accessories

Cover for hidraulic trolley

Wooden crate Package

