



Disclaimer:
Images are purely
indicative. The features
and the accessories of the
products may differ
according to model and
application.

Characteristics

External Dimensions (WxDxH)	cm 260x75x85
Height adjustable Feet	
Structure	Completely in stainless steel; work surface with a high, profiled surround and a seamlessly welded basin. Table carrying capacity: max. 200 kg. Basin cm 40 x 50 x 20 (W/D/H)
Weight	110 kg

Technical data

Package	Strong cardbox on wooden pallet <i>Dim. cm. 275 x 80 x 100; Kg.120</i> <i>Wooden crate on wooden pallet available</i>
----------------	---

Accessories

Hydroaspirator	
Body support	
Steel table for instruments support	
Rubber neck support	
Stainless steel neck support	
Knee-operated mixing tap	
Wooden crate Package	