



Disclaimer:  
Images are purely  
indicative. The features  
and the accessories of the  
products may differ  
according to model and  
application.

## Characteristics

<b>External Dimensions (WxDxH)</b>	260x80x85 cm
<b>Structure</b>	<p>Completely in stainless steel, the Autopsy table is mounted on a center support column. The work surface features a high, profiled surround and a seamlessly welded basin is affixed to the end of the table, with a negative inclination arranged diagonally to the drain. All power connections are built into the base. Included:</p> <ul style="list-style-type: none"> <li>• Central support column with one access doors</li> <li>• Basin W/D/H 40 x 50 x 20 cm</li> <li>• 2" Outlet with strainer</li> <li>• Knee-operated mixing tap far cold and warm water</li> <li>• 3 m shower hose with hand sprinkler</li> <li>• 2 splash proof electrical outlets</li> <li>• Electrical and pumping connections are preinstalled</li> </ul>
<b>Weight</b>	200 kg

## Technical data

<b>Package</b>	<p>Strong cardboard on wooden pallet  <i>Dim. cm. 276x101x144; Kg.267</i>  <i>Wooden crate on wooden pallet available</i></p>
----------------	---

## Accessories

<b>Hydroaspirator</b>
<b>Body support</b>
<b>Steel table for instruments support</b>
<b>Rubber neck support</b>
<b>Stainless steel neck support</b>
<b>Grinder</b>
<b>Wooden crate Package</b>