



Disclaimer:
Images are purely
indicative. The features
and the accessories of the
products may differ
according to model and
application.

Characteristics

External Dimensions (WxDxH)	cm 220x75x85
Rollers	Castors with brakes
Structure	Completely in stainless steel. The work surface features a high, profiled surround and a seamlessly welded basin is affixed to the end of the table, with a negative inclination arranged diagonally to the drain.
Weight	65 kg

Technical data

Package	Strong cardbox on wooden pallet <i>Wooden crate on wooden pallet available</i>
----------------	---

Accessories

Body support
Steel table for instruments support
Rubber neck support
Stainless steel neck support
Liquid bin
Wooden crate Package