

MORTUARY HOLDING TABLE

MOD.C

Cod. TVA007



Disclaimer:
Images are purely
indicative. The features
and the accessories of the
products may differ
according to model and
application.

Characteristics

External Dimensions (WxDxH)	cm 210x75x80
Structure	It is entirely made of AISI 304 stainless steel. This body preparation table features a top with a perimeter border to contain liquids. Equipped with drain pipe and a small stainless steel basin to collect liquids. Maximum load capacity: 125 kg

Accessories

Castors (H +50 mm)

Wooden crate Package