

AUTOPSY TABLE MOD.B VERTICALLY ADJUSTABLE AND 90° ROTATING

Cod. TVA011



Disclaimer:
Images are purely
indicative. The features
and the accessories of the
products may differ
according to model and
application.

Characteristics

External Dimensions (WxDxH)	260x85x85/110 cm
Structure	Made of AISI 304 stainless steel. With a reinforced and shaped structure, it consists of a platform with perimeter edging for the containment of liquids, without sharp edges and suitably inclined to allow easy flow of water towards the drain. The top includes a series of 5 removable perforated panels supporting the body of the deceased, the outflow of washing water and the suction of exhausted air directly from the surface. The provided stainless steel sink is equipped with taps and fittings with mixer for hot and cold water operated by elbow, wash hand shower extendable over the entire surface of the table. The table is equipped even with an irrigation system along with all the plane through water ejectors placed on the side opposite to the wastepipe.
Equipped with	· Full-floor irrigation system using water ejectors placed on the side opposite the drain. · Electric shredder-waste disposer. · Hydro-aspirator.
Supporting structure and table elevation system, 90° rotating version	The load-bearing element consists of a central structure made entirely of stainless steel adequately reinforced for weight load, equipped with removable panels to allow access to internal water and electrical connections. The platform is supported directly on an electric linear actuator that allows its elevation from a minimum of 850 mm to a maximum of 1100 mm. It is operated by means of a special switch with ascent and descent controls and a magnetothermal panel with a differential for protection. The platform is rotatable up to 90°.
Stainless steel basin dimensions	40x25x50 cm
Max. load capacity	200 kg

Technical data

Voltage	230V/1~/50hz
----------------	--------------

Accessories

Alimentation 220 V/1~/60 Hz
3-piece body support for autopsy table, 800 mm depth
Perforated organ tray for 800 mm worktop
Smooth instrument tray for 800 mm worktop
Perforated organ tray for 850 mm worktop
Smooth instrument tray for 850 mm worktop
Multi-position headrest
Wooden crate Package